

OXY-GUARD **HAND SANITIZER Gel**

ALCOHOL ANTISEPTIC 80% TOPICAL SOLUTION

Oxyguard hand sanitiser gel is a perfumed, quick drying hand sanitizing gel, highly effective against a wide range of Bacteria & Viruses. A combination of Hydrogen Peroxide and Ethyl alcohol, Oxyguard is pure and effective during high-risk environment. Hydrogen peroxide works as a flash disinfectant by destroying essential components of germ cells, and can deactivate a wide range of microorganisms, including bacteria, viruses, fungi, and spores. Ethyl alcohol exhibit rapid broad-spectrum antimicrobial activity against bacteria, viruses, and fungi and in amalgamation significantly increase product efficacy.



Salient features :

- Powerful germ killing efficiency. Kills 99.99% of germs
- Cleans & freshens hands.
- Effectively removes a wide range of stains.
- Moisturizing formulation prevents skin dryness.
- Effective against enveloped and non-enveloped viruses
- Enticing fragrance note.

Composition:

Ethyl Alcohol IP(Denatured with DEP): 80% v/v
 Glycerine IP (99%): 1.4% v/v
 Alkyl Acrylate cross polymer: 0.6% v/v
 Hydrogen peroxide IP 3%: 0.125% v/v
 Triethanolamine: 0.2% v/v
 Perfume: 0.3% v/v
 DM Water: qs to 100%

Specification :

Form : Liquid gel
 Colour : Transparent
 Odour : Perfumed
 Sp.Gravity : 0.88 – 0.89
 pH : 6.3 to 6.5
 Boiling point range : 80°C – 83°C
 Melting point : Not tested
 Flash point : 28°C
 Methyl alcohol content : Nil

Instructions for use :

- Squirt as needed into your palms to completely spread both hands
- Interlace your fingers with the palm against palm. Give special attention to finger tips, beneath the nails and thumbs.
- Rub hands together, both palms to palms and fingers to fingers.
- Rub both your fists
- Hold the left thumb with your right hand and rub repeatedly. Reciprocate likewise the right thumb
- Keep rubbing until your hands are dry.

Precautions :

Always read the label and use only as directed
 Always use undiluted
 Topical solution – strictly for external use
 Avoid contact with eyes. In case of eye contact flush with water
 Do not mix with soap or other chemicals
 Discontinue use if irritation occurs
 Do not use after the expiry date.

Item	volume	Delivery Form
Oxyguard	100 ml	40 units / carton
Oxyguard	500 ml	24 units / carton
Oxyguard	1000 ml	12 units / carton
Oxyguard	5000 ml	4 units / carton
Oxyguard	25000 ml	Discrete



Incon Labs

Kormala, Kuttichira P.O, Chalakudy, Thrissur- 680 724, Kerala State.

Customer Care No: +91 8113800999

Email: labs@incon.net.in

Mfg. Lic. No: 43/25/2020



Keep out of reach of Children



Flammable